PROTOCOLLI CHIRURGICI
SURGICAL PROCEDURES

B&B DENTAL
implant company

INDICE

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INTRODUCTION

This manual is designed for being used by clinicians who have undergone at least basic surgical and in-clinic implant training. All the information about the Duravit system show to dentists and specialists essential steps regarding implant planning and surgical procedures.

B.& B. Dental has an interest to keep doctors up-dated over the latest trends and treatment techniques of implants through a continued education. Function, beauty and biology in perfect harmony.

Our aim is to provide you a wide range of implant solutions. We develop products and solutions in order to make your job as simple as possible even in case of complex surgery.

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
DOUBLE THREAD BODY

- Thread of increased depth and highly sharp
- Easy insertion and bone condensing
- Higher primary stability

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BACK-TAPERED CORONAL DESIGN WITH MICRO RINGS

- Optimal soft tissue support
- Maximum alveolar bone volume
- Less crestal resorption

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APICAL BLADE

- Penetrate small diameter preparations
- Optimal anchorage

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











MICROGROOWING COLLAR

- Primary stability increased
- Easy insertion implant
- Reduced prosthetic load
- Easy gingival tissues healing

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PERFECT ANATOMICAL DESIGN LIKE A RADICULAR STRUCTURE

- Innovative macromorphology, designed for a primary high stability

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




DURAVIT
● ø3,5

PROCEDURA STANDARD DRILLING PROTOCOL - STANDARD

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

The Bone Tap is available as an optional step within the drilling sequence in case of dense bone (D1).

Codice colore	● ø 3,5	Color Code
Diametro Reale	ø 4	Real diameter
ref: TAP-35		

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






An efficient and atraumatic preparation of the implant site must be realized through a step-by-step procedure using drills of different diameters. All drilling phases of the bone tissue should be carried out under profuse external irrigation with saline solution or sterile distilled water and with an intermittent drilling technique to prevent overheating of the bone and to create a draining effect for an efficient removal of bone tissue.

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Cut and open the flap widely enough to have a complete view of the inclination and width of the bone ridge. Use the LANCE DRILL to pierce the cortical bone and prepare the implant head site, facilitating the insertion and the direction of the following $\varnothing$ 2.1 depth drill.

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Standard Metal stoppers

They are press-fitted in drills or compactors ($\varnothing$ 2.1/ $\varnothing$ 3/ $\varnothing$ 3.5/ $\varnothing$ 4/ $\varnothing$ 4.5), allowing an easy and quick perforation of the planned implant site.

Note: $\varnothing$ 5 drill and compactor-expander mount a different type of stopper.

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Depth drill $\varnothing$ 2.1 and $\varnothing$ 3

By using the proper stopper, drill paying attention only to the inclination.

Note: Parallel pins can be inserted into the implant site, facilitating the direction of the subsequent drilling.

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*Standard Drill $\varnothing 3.5$ and $\varnothing 4$
Keep drilling the implant site reaching the appropriate depth and width.*

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*IMPLANT EXTRACTION
Pull out the implant from the packaging using the specific driver manual or contra-angle, which can be both manual and for the universal contra-angle.*

Important Note: The driver shall be provided with the spring.

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*IMPLANT INSERTION
Screw the implant into the prepared site with the specific driver, until you feel good stability, then remove it.*

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
*FINAL IMPLANT INSERTION
Ratchet and implant driver will be used in order to complete the implants placing 1-1,5 mm inside the site.*

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Duravit bone Compactor-Expanders are a good alternative of osteotomes for the expansion and condensation of the bone in particular in the maxilla, during preparation of implant site. They are also an alternative to the Summers maxillary sinus lift technique. Compactor-Expanders increase implant success improving primary stability of the implant and maintenance of bone density, increasing fixation. They are driven into the bone manually with a straight surgical driver or with torque ratchet, decreasing the trauma typical of osteotomes.

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Cut and open the flap widely enough to have a complete view of the inclination and width of the bone ridge. Use the LANCE DRILL to pierce the cortical bone and prepare the implant head site, facilitating the insertion and the direction of the following $\varnothing 2.1$ depth drill.

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
STANDARD METAL STOPPERS

They are press-fitted in drills or compactors ($\varnothing 2.1/ \varnothing 3/ \varnothing 3.5/ \varnothing 4/ \varnothing 4.5$), allowing an easy and quick perforation of the planned implant site.

Note: $\varnothing 5$ drill and compactors-expander mount a different type of stopper.

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DEPTH DRILL $\varnothing 2$ AND COMPACTOR-EXPANDER $\varnothing 3$

By using the proper stopper, drill paying attentions only to the inclination. Then carry on using the Compactor-Expander $\varnothing 3$, connected to the manual straight key, compact and prepare the site.

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COMPACTOR-EXPANDER Ø 3.5 AND Ø 4

Keep on compacting the implant site until we feel a good resistance of the tools, in order to guarantee a great primary stability of the implant.

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IMPLANT EXTRACTION

Insert the ratchet mounter provided with the spring into the straight manual key and pull out the implant from the packaging.

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IMPLANT INSERTION

Screw the implant into the prepared site with the ratchet mounter connected to the straight manual key and remove it when you feel good stability.

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FINAL IMPLANT INSERTION - OPZIONE 1

Ratchet and implant driver will be used in order to complete the implants placing into the site.

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FINAL IMPLANT INSERTION - OPZIONE 2

Straight manual key and implant driver will be used in order to complete the implants placing into the site.

PROCEDURA SEPOLTA E TRANSMUCOSA

ONE AND TWO-STAGE PROCEDURES

PROCEDURA SEPOLTA TWO-STAGE PROCEDURE

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Connect the hex screwdriver with the cover screw, placed on the top of the plastic plug

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Using light finger force (5 Ncm) seat it inside the implant.

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Suture the flaps carefully and tighten them together.

PROCEDURA TRANSMUCOSA ONE-STAGE PROCEDURE

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Connect the screwdriver with the transmucosal screw placed on the bottom of the plastic plug.

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Suture carefully the soft tissue flaps around the transmucosal screw.

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UNDER CRESTAL BONE

It is advisable to place always the implant minimum 1 mm under crestal bone:

- the healing of the bone tissue is facilitated, reducing the risk of bone loss;
- soft tissues are maintained high in order to guarantee an optimal esthetic effect for long time.

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BONE LEVEL

It is necessary to place the implant at the bone level in the following conditions:

- there is enough gingiva for closing the platform switching;
- whenever it is inserted in mandible and the nerve is very close to it.

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NEVER OVER THE BONE LEVEL

Never place the implant over the bone level because it could occur the reduction of the soft tissues, preventing the use of abutment with platform switching.

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
IMPLANT TREATMENT FEASIBILITY

this study takes different elements into consideration:


- A questionnaire for patients to reveal eventual health problems/medications which could have an impact on the treatment success, alcohol, tobacco or drugs consumption, general dental hygiene...
- An oral examination which will give details about the mouth opening, the line of the patient's smile (is it a gingival smile ?) the restorative height and the available bone volume and the occlusion type...
- Biological tests (Glycemy...)
- A complete X-Ray file showing the available bone volume.
- Full study models with the two dental arches in occlusion.
- An implant treatment procedure cannot be started without a complete cleaning of all the patient's infectious cases.
- Fill in the informed consent with patient.

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DENTAL PANORAMIC RADIOGRAPH

- This X-ray plate shows an enlargement of 1:25 or 1:3 depending on the device being used.
- Thanks to the template supplied with the surgical kit the doctors can identify easily implant lengths, helped by implant profile increased of 25%.
- It enables to study the bone height in the edentulous zone but it does not indicate the crest thickness.
- This examination is essential but not sufficient for a pre-implant evaluation.

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
Arcata intera
Full Arch

C.T. DENTALSCAN - CONE BEAM


- *It is the most performant and the most reliable technique in the preimplant evaluation.*
- *It enables a precise anatomical study of the edentulous area without any distortion.*
- *It points out anatomical obstacles (Nasal sockets, sinus, mandibular canal)*
- *It enables a real size morphological study of the toothless area: available bone height, labio-lingual thickness, jaws obliquity.*
- *An analysis of the bone structure is obtained: condition of the spongy tissue and thickness of the cortical bone*
- *The documents are in real scale (ratio 1/1), so it is possible to superimpose the implant template on real scale section.*

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Care must be taken to avoid the inferior alveolar nerve and the mental foramina in the premolar region, since the mandibular nerve is often inclined coronally in this area.

The location of the maxillary sinus and nasal floor must be positively identified to avoid the inadvertent penetration with a reamer or an implant.

In general, at least 2.0 mm of bone should separate the apex of the implant from the mandibular canal.

Care must be taken to avoid the penetration of the submandibular fossa which is located below the mylohyoid line and particularly the sublingual space in the anterior mandible where the sublingual artery is located.

Inadvertent penetration of these lingual plates may be avoided by directing appropriately the pilot bur and reamer burs toward the buccal and by monitoring the area with digital contact while drilling.

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

IMPLANT SIZE SELECTION

- *The appropriate implant length and width depends upon the available bone and the expected occlusal loads.*
- *In general, choose the widest and the longest implant possible.*
- *Panoramic and periapical radiographs as well as diagnostic models and a clinical examination are used to determine if enough mesio-distal space and vertical bone height exist to place a Duravit implant safely and appropriately in a proposed site.*
- *A transparent template, which depict implant outlines of actual size and 25% of actual size, is helpful in selecting an appropriate implant. Since radiographs are not necessarily precise representations, knowledge of their magnification must be considered while using them to determine an appropriate implant choice.*

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
The chart contains indicative values only.

Actual clinical conditions and the clinician's assessment of the patient should be the main criteria for choosing the size of an implant for a specific area.

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The surgical organizer is used for the secure storage and sterilization of the surgical instruments.

It's made of a highly shock-proof thermoplastic, suitable for frequent sterilization in the autoclave. Autoclaving at a temperature not over 121 °C/250 °F is recommended.

KIT COMPLETO
COMPLETE KIT

KIT CHIRURGICO WIDE WIDE SURGICAL KIT

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OVERVIEW OF SURGICAL INSTRUMENTS

Paralleling pin is inserted into the prepared implant site after using drill $\varnothing$ 2.1 or 3, facilitating the direction of the subsequent drilling.

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Tissue punch allows to pierce the mucosa according to the selected implant diameter, if it is used with a contra angle set at a low speed.

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Extender drill increases operating length of the drills during the surgery.

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Lance drill marks out and creates the insertion point penetrating the cortical bone.

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WIDE DRILLS available in 4 diameters ($\varnothing$ 5.5/ $\varnothing$ 6/ $\varnothing$ 6.5/ $\varnothing$ 7) prepare the implant site. Useful for immediate placement in molar extraction socket.

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Duravit bone **COMPACTOR-EXPANDER** are an alternative to osteotomes for the expansion and condensation of the atrophic mandible and maxilla in preparation for dental implant insertions. Compactorexpander are also an alternative to the Summers maxillary sinus lift technique.

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COUNTERSINK is used to create a shoulder area for the implant. Indicated especially in hard bone conditions to avoid compression around the cortical bone.

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BONE TAPS may be used to prepare the implant thread profile into the implant site reducing the bone pressure.

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The drill **STOPPERS** ensure to the surgeon simple and accurate depth control, quick and easy assembling as well. The laser mark number identifies the length.

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RATCHET MOUNTERS are connected to the ratchet are used to place the implants "3P" and "WIDE". It's available LONG and SHORT.

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
OVERVIEW OF SURGICAL INSTRUMENTS

CONTRA-ANGLE implant **DRIVER** connected to the contra-angle is used to the initial placement of 3P and WIDE implants; reached the optimal resistance, screw with the ratched mounters.

RATCHET MOUNTERS $\varnothing 3$ are connected to the ratchet and used to place the implants " $\varnothing 3$ ". Available long and short.

HANDPIECE DRIVER 1.25 Hexagon connected to the contra-angle is used to set all the components screws such as Cover screw, Healing screws, Prosthetic screws.

MANUAL SCREW DRIVERS 1.25 Hexagon are used to screw all the prosthetic components such as Cover screw, Healing screws, Prosthetic screws. Available long and short.

TORQUE SCREW DRIVERS 1.25 Hexagon, connected to the torque ratchet, is used to set with predefined torque force, all the components screws such as cover screw, healing screws, prosthetic abutments. Available long and short.

TORQUE RATCHET allows the clinician to apply accurately the recommended torque for surgery and prosthetics, adjustable between 0 and 50 nc2.

FIX RATCHET allows the clinician to place the implants inside the prepared implant site, without using any predefined torque force.

STRAIGHT MANUAL KEY allows the clinician to prepare the implant site by using the compactor-expander and then to insert the implant. It is especially used in the maxillary.

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CLEANING AND MAINTENANCE OF THE TORQUE RATCHET

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For ensuring that the torque ratchet works properly, it must be always taken apart and cleaned immediately after each surgery. (See picture above) Otherwise, if blood penetrates into the instrument, it will stick the ratchet bolt into the ratchet body, which means corrosion!

Disinfect, clean, and sterilize as described on the following pages. Please treat each component separately!

After reassembling the torque ratchet, its correct working must be checked, using the implant driver. Holding the driver tightly, the ratchet should turn. The same control must then be carried out for the opposite direction.

PLEASE NOTE:

We recommend to keep always at disposal another sterilized torque ratchet during the surgery.

Check promptly that instruments work in a correct way before the surgery!

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CLEANING PHASES

Instruments with high cutting performance are a basic requirement for successful implantation. Therefore the following should be remembered:

- *Cleansing.*
- *Dirty instruments should be placed in a bowl with a cleaning solution in order to disinfect them, avoiding that the residues of blood coagulate.*
- *With a brush remove immediately from the instruments residues of blood, saliva, tissue after surgery. If residues adhere to the instruments and dry on them, they may lead to corrosion.*
- *Disinfection.*
- *Clean them accurately using an ultrasonic device and a disinfectant liquid to avoid the corrosion of the drills. Do not use chlorine-based or peracetate disinfectants.*
- *An accurate rinsing preferably with distilled water and a subsequent drying conclude the preparation phase for steam sterilization.*

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STERILIZATION PHASES

After cleaning, all surgical instruments must be sterilized:

- Only when the instruments are cleaned, disinfected and dried, they can be sorted into the surgical organizer, which has to be wrapped in sterilizable envelopes.

NOTE: The sterilizable envelopes are easily commercially available.

- For a long duration of the surgical organizer sterilization must be preferably performed in steam autoclave at a temperature of 121 °C/ 250 F°.

IMPORTANT:

- Chemical sterilization is not recommended since this procedure can damage the plastic surface.
- No hot air sterilization since the high temperature (appr. 180 °C) would make the plastic cassette melt.
- No sterilization in spherical sterilizers (temperature up to 300 °C).
- In order to avoid damages to the surgical cassette during steam sterilization, the surgical cassette has to be placed correctly in the autoclave.

Via San Benedetto, 1837 - 40018 San Pietro in Casale (BO) Italy
Tel. +39 (0) 51.81.13.75 - Fax +39 (0) 51.666.94.00
info@bebdental.it - www.bebdental.it

Certified quality system
UNI EN ISO 13485