

# HEALING COMPONENTS

## HEALING SCREW (Titanium grade 5)

Provided in the implant package. Use when you want to cover completely the implant after the insertion.  
The implant is reopened 3-6 months, followed by the use of healing screws.



016013



### TIGHTENING:

Insert the healing screw into the implant and tighten with only light finger force.

## HEALING SCREWS (Titanium grade 5)

Each package contains 1 piece.

These components are used to rehabilitate soft tissue on the implant in order to insert the final prosthetic abutment later on.

They are available in different heights and in 2 different diameters profile:

- Ø 5 indicated for anterior area
- Ø 6 indicated for posterior area



INN-6010

INN-6011

INN-6012

INN-6014

INN-6015



INN-6020

INN-6021

INN-6022

INN-6023

INN-6024



### TIGHTENING:

Insert the healing screw into the implant and tighten with only light finger force.

## HEALING SCREWS (TITANIUM GRADE 5)

Each package contains 1 piece.

These components are used to rehabilitate soft tissue around the implant in order to insert the final prosthetic abutment later on.

